



Recombinant Human immunodeficiency virus type 1 group M subtype B Gag polyprotein (gag)

Product Code	CSB-MP758482HDAK
Abbreviation	gag
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q73367
Product Type	Recombinant Protein
Immunogen Species	Human immunodeficiency virus type 1 group M subtype B (strain 89.6) (HIV-1)
Purity	>85% (SDS-PAGE)
Sequence	GARASVLSG GELDRWEKIR LRPGGKKKYK LKHIVWASRE LERFAVNPSL LETSEGCRQI LGQLQSSLQT GSEELKSLYN TVATLYCVHQ RIEVKDTKEA LDKIEEEQNK SKKKAQAAAA DTGNSSQVSQ NYPIVQNIQG QMVHQAISPR TLNAWVKVVE EKAFSPEVIP MFSALSEGAT PQDLNTMLNT VGGHQAAMQM LKETINEEAA EWDRLHPVQA GPVAPGQMRE PRGSDIAGTT STLQEIQIGWM TNNPIPVGE IYKRWILGL NKIVRMYSPTS SILDIKQGPK EPFRDYVDRF YKTLRAEQAS QEVKNWMTET LLVQANPDC KTILKALGPG ATLEEMMTAC QGVGGPGHKA RVLAEAMSQV TNSATIMMQR GNFRNQRKTV KCFNCGKEGH IAKNCRAPRK KGCWKCGKEG HQMKDCTERQ ANFLGKIWPS HKGRPGNFLQ SRPEPTAPPE ESRFRGEETT TPSQKQEPID KELYPLASLR SLFGNDPSSQ
Source	Mammalian cell
Target Names	gag
Protein Names	Recommended name: Gag polyprotein Alternative name(s): Pr55Gag Cleaved into the following 6 chains: 1. Matrix protein p17 Short name= 2. MA 3. Capsid protein p24 Short name= 4. CA 5. Spacer peptide p2 6. Nucleocapsid
Expression Region	2-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.