



# Recombinant *Ashbya gossypii* Actin-related protein 4 (ARP4)

<b>Product Code</b>	CSB-MP758714DOT
<b>Abbreviation</b>	ARP4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q754G5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSNSALQVYG GDEITAVVID PGSFTTNIGY SGTDCPQAIL PSCYGKYTEG EKDELFFSEQS IGLPRKDYEI HNIVQNGEVV DWEKAEKQWD WAIRSELRFE TNSGMPALLT EPIWNSEENR KKSLEVLLES MDFSACYLVP TATAVSFAMG RPTCLVVDIG HDVTSVCPVV DGMTLSKSSM RSYIAGSLLN ELIRSQLAPR KVIPLFQVAQ RRPVFMERKF DYEIHPSLQK FVNERQFFQE FKETMLQVAP TSISKFKSEI ETTSKRSIEA PWGEELVYDS LQRLEFAEQL FTPDLSQFPE DWPISKDGVV ETWHNDYVPL KRNKPGTNVK DKEGTLDATP VPDENSVTSA DQPNNDNGKRN LEETTPDQKN EVSGLADLIY SSIMSTDVDL RTTLSHNVVI TGGTSSLPGL MDRISAELNR SLPALKFRML TSGQLRERQY QGWLGGSSILA SLGTFHQLWV GKQEYAEVGA DRLLKDRFR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ARP4
<b>Protein Names</b>	Recommended name: Actin-related protein 4 Alternative name(s): Actin-like protein ARP4 Short name= Actin-like protein 4
<b>Expression Region</b>	1-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.