



Recombinant *Ashbya gossypii* 40S ribosomal protein S1 (RPS1)

Product Code	CSB-EP758731DOT
Abbreviation	RPS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q756Y3
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	AVGKNKRLS KGKKGLKKKV VDPFTRKEWY DIKAPSTFEN RNVGKTLVNK STGLKNAADF LKGRVVEVCL ADLQGSEDHS FRKVKLRVDE VQGKNLLTNF HGLDFTTDKL RSMVRKWQTL IEANVTVKTA DDYVLRVFAI AFTRRQSNQI KKTSYAQSSH IRAIRKWISE ILTREVQNAT LAQLTSKLIP EVINKEIENA TKDIFPLQNV HIRKVKLLKQ PKFDL GALMT LHGESSGEEK GKKVAGGFKD EILETV
Source	<i>E.coli</i>
Target Names	RPS1
Protein Names	Recommended name: 40S ribosomal protein S1
Expression Region	2-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.