



Recombinant *Ashbya gossypii* Coenzyme Q-binding protein COQ10, mitochondrial (COQ10)

Product Code	CSB-YP758771DOT
Abbreviation	COQ10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75CC1
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	SRRTFL GFTGGDTKEQ RYILKRVFNA PLHYVYPAVS EVSLYKLFIP YCTDSFVNKR RPGDNMPTEA GLRVGFQQYD ETFVCRVDCT TLPGNQRSVV AESLAHHLFE TLHTQWLLSP HPTRPDASVV ELILRFKFKS QLYNSVSSIF GTRVTQVVMK AFEKRVFQLR KEAMDRPPTG DAGH
Source	Yeast
Target Names	COQ10
Protein Names	Recommended name: Coenzyme Q-binding protein COQ10, mitochondrial
Expression Region	25-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.