



Recombinant *Ashbya gossypii* Protein SWT21 (SWT21)

Product Code	CSB-EP758785DOT
Abbreviation	SWT21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75E79
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MNILAHTGDI FHGRGKLHAS LKEQQVWDDL RRPGHDVPQL PRPQGYTASQ PIYKRLPVVC RGLTWSFDGS SLVAVHEDCG IRQYAVCEEQ GPVMAPVQRF FKPQAIISHA VHPQNSLFQA DSSLNVVVVS CKDMPYQLYP LHSNPDQPAR ALQTYNSASE STEEFRAAYA MDFLDDAHFL AGSTRNTVSL FDASRRESIW VSSATRRGCG RGQHKAIIVSC FDEANERRDC VRYAGTYRSE LLRIDPRTSS ISSWKAIQDK SRGIYQVLQS NNGHYFYIMK RYSNSITVLD SRKSMAIVNQ FKLPFQTHAQ KLRGTYSTRH GLLIGDEAGR ILQWAPEVVE FGGIDRSGSV PSYGIGETGI THVCAEGSRI NIIKPNADS DIFAISYSPD KYGEDVSCNA GILLMRL
Source	E.coli
Target Names	SWT21
Protein Names	Recommended name: Protein SWT21
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.