



Recombinant *Ashbya gossypii* Mediator of RNA polymerase II transcription subunit 20 (SRB2)

Product Code	CSB-MP758791DOT
Abbreviation	SRB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75F10
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MKRTAVIFIE KATPSTITQF HDILSTHVLA IQEKWSFELK TFRSSIKNLP PSDTKVLVSL QLTHRDNQTV VIKNQSAIVT GQHSTDALTS NGCSSGFPEP FDNILTSLKLS NIWTQRQSTK GNFGTTYKTS ELIIRASNVF SSSGFKGLLL EIECTDPVSA EEFDRRVANI RAMLSEIDIN DYKLNKDEM N EGKPVFLCDL AYQYVKVLD
Source	Mammalian cell
Target Names	SRB2
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 20 Alternative name(s): Mediator complex subunit 20
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.