



# Recombinant Human MTSS1-like protein (MTSS1L)

<b>Product Code</b>	CSB-EP758826HU
<b>Abbreviation</b>	MTSS1L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q765P7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	METAEKECGA LGGLFQAIVN DMKSSYPIWE DFNSKATKLH SQLRRTTVLAA VAFLDAFQKV ADMATNTRGA TRDIGSALTR MCMRHRSIET KLRQFTNALL ESLINPLQER IEDWKKAANQ LDKDHAKEYK RARHEIKKKS SDTLKLQKKA RKELLGKGD LQPQLDSALQD VNDMYLLLEE TEKQAVRRAL IEERGRFCTF ITFLQPVVNG ELTMLGEITH LQGIIDDLVV LTAEPHKLPP ASEQVIKDLK GSDYSWSYQT PPSSPSSSSS RKSSMCSAPS SSSSAKGGGA PWPGGAQTYS PSSTCRYRSL AQPATTTARL SSVSSHDSGF VSQDATYSKP PSPMPSDITS QKSSSSASSE ASETCQSVSE CSSPTSDWSK VGSHEQPSGA TLQRRKDRVE LLRDTEPGPA SGGTLGPSGE EAPRPRMSPA TIAAKHGEEV SPAASDLAMV LTRGLSLEHQ KSSRDSLQYS SGYSTQTTTP SCSEDTIPSQ GSDYDCYSVN GDADSEGPPE FDKSSTIPRN SNI AQNYRRL IQTKRPASTA GLPTAGLPTA TGLPSGAPPV VATIRRTPST KPTVRRALSS AGPIPIRPPI VPVKTPTVPD SPGYMGPTRA GSEECVFYTD ETASPLAPDL AKASPKRSL PNTAWGSPSP EAAGYPGAGA EDEQQQLAAN RHSLVEKLGE LVAGAHALGE GQFPFPTALS ATPTEETPTP PPAATSDPPA EDMLVAIRRG VRLRRTVTND RSAPRIL
<b>Source</b>	E.coli
<b>Target Names</b>	MTSS1L
<b>Protein Names</b>	Recommended name: MTSS1-like protein Alternative name(s): Actin-bundling with BAIAP2 homology protein 1 Short name= ABBA-1
<b>Expression Region</b>	1-747
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.