



Recombinant Human MTSS1-like protein (MTSS1L)

Product Code	CSB-MP758826HU
Abbreviation	MTSS1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q765P7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	METAEKECGA LGGLFQAIVN DMKSSYPIWE DFNSKATKLH SQLRTTVLAA VAFLDAFQKV ADMATNTRGA TRDIGSALTR MCMRHRSIET KLRQFTNALL ESLINPLQER IEDWKKAANQ LDKDHAKEYK RARHEIKKKS SDTLKLQKKA RKELLGKGD L QPQLDSALQD VNDMYLLLEE TEKQAVRRAL IEERGRFCTF ITFLQPVVNG ELTMLGEITH LQGIIDDLVV LTAEPHKLPP ASEQVIKDLK GSDYSWSYQT PPSSPSSSSS RKSSMCSAPS SSSSAKGGGA PWPGGAQTYS PSSTCRYRSL AQPATTTARL SSVSSHDSGF VSQDATYSKP PSPMPSDITS QKSSSSASSE ASETCQSVSE CSSPTSDWSK VGSHEQPSGA TLQRRKDRVE LLRDTEPGPA SGGTLGPSGE EAPRPRMSPA TIAAKHGEEV SPAASDLAMV LTRGLSLEHQ KSSRDSLQYS SGYSTQTTTP SCSEDTIPSQ GSDYDCYSVN GDADSEGPPE FDKSSTIPRN SNIAQNYRRL IQTKRPASTA GLPTAGLPTA TGLPSGAPPG VATIR RTPST KPTVRRALSS AGPIPIRPPI VPVKTPTPVD SPGYMGPTRA GSEECVFYTD ETASPLAPDL AKASPKRLSL PNTAWGSPSP EAAGYPGAGA EDEQQQLAAN RHSLVEKLGE LVAGAHALGE GQFPFPTALS ATPTEETPTP PPAATSDPPA EDMLVAIRRG VRLRRTVTND RSAPRIL
Source	Mammalian cell
Target Names	MTSS1L
Protein Names	Recommended name: MTSS1-like protein Alternative name(s): Actin-bundling with BAIAP2 homology protein 1 Short name= ABBA-1
Expression Region	1-747
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.