



Recombinant Varicella-zoster virus Transcriptional transactivator IE63 (ORF63)

Product Code	CSB-EP758863VAQ-B
Abbreviation	ORF63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q77NN7
Product Type	Recombinant Protein
Immunogen Species	Varicella-zoster virus (strain Oka vaccine) (HHV-3) (Human herpesvirus 3)
Purity	≥85% (SDS-PAGE)
Sequence	MFCTSPATRG DSSESKPGAS VDVNGKMEYG SAPGPLNGRD TSRGPAGFCT PGWEIHPARL VEDINRVFLC IAQSSGRVTR DSRRLRRICL DFYLMGRTRQ RPTLACWEEL LQLQPTQTQC LRATLMEVSH RPPRGEDGFI EAPNVPLHRS ALECDVSDDG GEDDSDDDGS TPSDVIEFRD SDAESSDGED FIVEEESSES TDSCEPDGVP GDCYRDGDGC NTPSPKRPQR AIERYAGAET AEYTAALKALT ALGEGGVDWK RRRHEAPRRH DIPPPHGV
Source	E.coli
Target Names	ORF63;
Protein Names	Recommended name: Transcriptional transactivator IE63 Alternative name(s): Immediate-early protein 63 Short name= IE63
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.