



Recombinant Staphylococcus aureus Endoribonuclease MazF (mazF)

Product Code	CSB-MP758890SUL
Abbreviation	mazF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7A0D7
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain MW2)
Purity	≥85% (SDS-PAGE)
Sequence	MIRRGDVYLA DLSPVQGSEQ GGVRPVVIIQ NDTGNKYSPT VIVAAITGRI NKAKIPTHVE IEKKKYKLDK DSVILLEQIR TLDKKRLKEK LTYLSDDKMK EVDNALMISL GLNAVAHQKN
Source	Mammalian cell
Target Names	mazF
Protein Names	Recommended name: mRNA interferase MazF EC= 3.1.-.-Alternative name(s): Endoribonuclease MazF Toxin MazF
Expression Region	1-120
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.