



Recombinant Putative uncharacterized protein ydet (ydeT)

Product Code	CSB-YP759010EOD
Abbreviation	ydeT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7DBJ9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MSGYTVKPPT GDSNEQTQFI DYFNLFYSKR GQEQISISQQ LGNYGATFFS ASRQSYWNTS RSDQQISFGL NVPFGDITTS LNYSYSNNIW QNDRDHLLAF TLNVPFSHWM RTDSQSAFRN SNASYSMSND LKGGMTNLSG VYGTLLPDNN LNYSVQVGNT HGGNTSSGTS GYSTLNYRGA YGNTNVGYSR SGDSSQIYYG MSGGIIAHAD GITFGQPLGD TMVLVKAPGA DNVKIENQTG IHTDWRGYAI LPFATEYREN RVALNANSLA DNVELDETVV TVIPHTGAIA RATFNAQIGG KVLMTLKYGN KSVPFGAIVT HGENKNGSIV AENGQVYLTG LPQSGKLQVS WGNDKNSNCI VDYKLPEVSP GTLLNQQTAI CR
Source	Yeast
Target Names	ydeT
Protein Names	Recommended name: Putative uncharacterized protein ydet
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.