



Recombinant Drosophila melanogaster Zinc finger CCCH domain-containing protein 15 homolog (CG8635)

Product Code	CSB-MP759057DLU
Abbreviation	CG8635
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7JWR9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPPKKAPPGP SKKTEQKKKE KVIEDKTFGL KNKKGNKQQK FIQQVQKQVQ SGGQHPRQDG DKRKEEKEKK LAEQREMAI FKPVQTQKVE KGTDPKSVVC AFFKQGTCTK GDKCKFSHDL SQENKVEKRS IYVDMRDEDD PMTNWDDAKL KEVVEKKSSG EKQRPTDII CKFFLEAVEK SKYGWFWECP NGGKCIYRHA LPAGYVLKRD KKKEEKPTI SLVDLIEKER AALGPNQTRV TLESFLAWKK RKIAEKKAKL AAEERKKSD FSKGKQFGIS GREMFSFNPD LVDDGPMEEG DAAFDVYNRE DDDDNAVEFK ELDLAALSLA AKEVDGSGTI ASTNRLLDQA TEAAKTA AE DAADEDGPS SSAPANDAAP INKDLFV DLA GELDDLDLDD EDDD
Source	Mammalian cell
Target Names	CG8635
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 15 homolog Short name= DRG family regulatory protein 1 homolog
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.