



Recombinant *Drosophila melanogaster* Coiled-coil domain-containing protein 130 homolog (CG15084)

Product Code	CSB-MP759058DLU
Abbreviation	CG15084
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7K0F0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGERKGQNKY YPPDYDPKKG GLNKFQGTHA LRERARKIHL GIIIRFEMP YNIWCDGCKN HIGMGVRYNA EKTKVGMYYT TPVFKFRMKC HLCDNHFEIQ TDPGNLDYVI LSGARRQENR WDPLQNEQVV PETKEVQKRL FDDAMYKLEH QAKDAKAGAD ARPVLQKLVE RNMSVWDDSY MANSRLRAEF RQQKKEINGQ QELDRQLLAK SSLDIALLEPE TTQDREMAAL MKLQTKSALE RESEQRLELL MRPALPGATV TTFGGLKRQK VLNTQLQVQD LGIRRKLEE TTSSATNEKP ISLVGDYSSS DNDSNG
Source	Mammalian cell
Target Names	CG15084
Protein Names	Recommended name: Coiled-coil domain-containing protein 130 homolog
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.