



Recombinant Human Protein FAM194A

(FAM194A)

Product Code	CSB-BP759168HU
Abbreviation	FAM194A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L0X2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAHLRSPSGF GDPGKKDQKE SEEELEEEEE EEEVEEEEEEE VEEEEEEVEE EEEEVVEEEL VGEEQLEAP ETFSEEYLWK VTDIGDYDDD FPDVRPRLAS IVSPSLTSTF VPSQSATSTE TPSASPPSST SSHKSFPKIF QTFRKDMSEM SIDRNIHRNL SPGIPVSVQT EESWLQDLSD KVQSRKKASK EKAPEPECLAS KLREKWINP EESKLNILYE LEFKEDFITL FEPSLRTLPS IGPPSILAYK EESSNLGINF KDEEEETSPK CEFCGSDLRA FFSNVDVSSE PKGHASCCIA FQNLIDYIYE EQIKTKPPKA ELIAIDPHAA HGSEVDRLKA KEKALQRKQE QRMARHFAII SREQTHFSED DSKRLKTISY QLSVDIPEKQ IIDDIVDFDQ LRNSNMSIIC CDSRIACGKV VRNELLEKHY KHGSKFLTSTF PDGTTQIFYP SGNLAIIRVP NKVNGFTCIV QEDMPTNPAI LAVLDSSGRS SCYHPNGNVW VYINILGGQY SDQAGNRIRA WNWSNSITSS PFVSFKPVFL ALNRYIGVRI LEQDKISITF LAMGQQARIS VGTKVKLPNP EEIPILRYVS GDDLILLASL IKIRRLFHKL EGCNVFPSSQ VWEKCLKQPSY LSSLSLKLIA LCHSSGIKQD IMKTIRNIIN EEI
Source	Baculovirus
Target Names	ERICH6
Protein Names	Recommended name: Protein FAM194A
Expression Region	1-663
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.