



Recombinant Human Vacuolar fusion protein MON1 homolog B (MON1B)

Product Code	CSB-MP759169HU
Abbreviation	MON1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L1V2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEVGGDTAAP APGGAEDLED TQFPSEEARE GGGVHAVPPD PEDEGLEETG SKDKDQPPSP SPPPQSEALS STSRLWSPAA PENSPTCSPE SSSGGQGGDP SDEEWSRQRK HVFVLSEAGK PIYSRYGSVE ALSATMGVMT ALVSFVQSAG DAIRAIYAED HKLVFLQQGP LLLVAMSRTS QSAAQLRGEL LAVHAQIVST LTRASVARIF AHKQNYDLRR LLAGSERTLD RL LDSMEQDP GALLLGAVRC VPLARPLRDA LGALLRRCTA PGLALSVLAV GGRLITAAQE RNVLAECRLD PADLQLLLDW VGAPAF AAGE AWAPVCLPRF NPDGFFYAYV ARLDAMPVCL LLLGTQREAF HAMAACRRLV EDGMHALGAM RALGEAASF S NASSASAPAY SVQAVGAPGL RHFLYKPLDI PDHHRQLPQF TSPELEAPYS REEERQRLSD LYHRLHARLH STSRPLRLIY HVAEKETLLA WVTSKFELYT CLSPLVTKAG AILVVTKLLR WVKKEEDRLF IRYPPKYSTP PATSTDQAAH NGLFTGL
Source	Mammalian cell
Target Names	MON1B
Protein Names	Recommended name: Vacuolar fusion protein MON1 homolog B Alternative name(s): HSV-1 stimulation-related gene 1 protein HSV-I stimulating-related protein
Expression Region	1-547
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.