



Recombinant Human Hsp90 co-chaperone Cdc37-like 1 (CDC37L1)

Product Code	CSB-YP759170HU
Abbreviation	CDC37L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L3B6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEQPWPPPGP WSLPRAEGEA EEESDFDVFP SSPRCPQLPG GGAQMYSHGI ELACQKQKEF VKSSVACKWN LAEAQQKLGS LALHNSESLD QEHAKAQTAV SELRQREEEW RQKEEALVQR EKMCLWSTDA ISKDFVFNKSF INQDKRKDTE DEDKSESFMQ KYEQKIRHFG MLSRWDDSQR FLSDHPYLVC EETAKYLILW CFHLEAEKKG ALMEQIAHQA VVMQFIMEMA KNCNVDPKRC FRLFFQKAKA EEEGYFEAFK NELEAFKSRV RLYSQQSQSFQ PMTVQNHVPH SGVGSIGLLE SLPQNPDYLQ YSISTALCSL NSVVHKEDDE PKMMDTV
Source	Yeast
Target Names	CDC37L1
Protein Names	Recommended name: Hsp90 co-chaperone Cdc37-like 1 Alternative name(s): Hsp90-associating relative of Cdc37
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.