



# Recombinant Mouse G2/mitotic-specific cyclin-B1 (Ccnb1)

<b>Product Code</b>	CSB-BP004806MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P24860
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALRVTRNTK INAENKAKVS MAGAKRVPVT VTAASKPGLR PRTALGDIGN KVSEELQARV PLKREAKTLG TGKGTVKALP KPVEKVPVCE PEVELAEPEP EPELEHVREE KLSPEPILVD NPSPSPMETS GCAPAEYYLC QAFSDVILAV SDVDADDGAD PNLCSYVVD IYAYLRQLEE EQSVRPKYLQ GREVTGNMRA ILIDWLIQVQ MKFRLLQETM YMTVSIIDRF MQNSCVPKKM LQLVGVTAMF IASKYEEMYP PEIGDFAFVT NNTYTKHQIR QMEMKILRVL NFSLGRPLPL HFLRRASKVG EVDVEQHTLA KYLMELSMMLD YDMVHFAPSQ IAAGAFCLAL KILDNGEWTP TLQHYLSYSE DSSLPVMQHL AKNVVMVNCG LTKHMTVKNK YAASKHAKIS TLAQLNCTLV QNLSKAVTKA
<b>Source</b>	Baculovirus
<b>Target Names</b>	Ccnb1
<b>Protein Names</b>	Recommended name: G2/mitotic-specific cyclin-B1
<b>Expression Region</b>	1-430
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found, a constitutively expressed transcript and a cell cycle-regulated transcript, that is expressed predominantly during G2/M phase. The different transcripts result from the use of alternate transcription initiation sites.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.