



Recombinant Danio rerio Putative deoxyribonuclease TATDN1 (tatdn1)

Product Code	CSB-EP760660DIL
Abbreviation	tatdn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GML7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTNFRFIDIG INLTDPMFRG VYRGTQKHED DFAEVVERAL QVGVQKFIIT GGNLEDSRAA LTLHTREQF FSTVGCHPTR CSEFDDQGSQ QYLSSLLDLT VSNTQKVAV GECGLDFDRL EFCPKETQLR YFQLQFDLAE ASGLPMFLHC RNAHTEFIDI MRRNRQRQCVG GVVHSFDGSQ QDAAALLDLT LYIGINGCSL KTAENLEVMK SIPSDRLMIE TDAPWCGIKN THAGAKLIKT SFPTKKKWET GHCVKDRNEP CHIIQVLEVM AAVREEDPLD LAETIFNNTD TLFFKNRS
Source	E.coli
Target Names	tatdn1
Protein Names	Recommended name: Putative deoxyribonuclease TATDN1 EC= 3.1.21.-
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.