



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 13 (Os02g0161200, LOC_Os02g06584)

Product Code	CSB-MP760705OFG
Abbreviation	Os02g0161200, LOC_Os02g06584
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6H7U2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MLGPPRRGPA YKTKLCALWQ RGNCNRDTCs FAHGHDIRr PPSSRGAFTH HPGRRDYRAG DFRGRIDRRF SPRRRHSPGR ESRGRPLYD RRPSSRERDS SYSRSPSRKS ERRHEKKTDD GETNSSRSLs LSDNNDEKKK DKFSSGDEKE DHEKQLKQIR LDMEALRDDK TQMEVILDEK IDEVRKISSK VNDLEVQLRR EKDECHRMTS KMKKFIKAHA RFLKAQEEVK RSQARFERLG DLLASDILKR GANEEGSSVN EDLNERSPNT AATKKRSIPY STSEAKAVK KRRERDSDTM TRSDKYRSDV TDFDKTSKGT EATKSLYLKK KLWEDEKSKL GANIFTEKVK GSPVRHVLPS TGMAAHAIDD LNEAIELEDR HESIDALLEN DADDKTRSPA IPLQPPPVVQ NAYEQYEGDD EEVDVE
Source	Mammalian cell
Target Names	Os02g0161200
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 13 Short name= OsC3H13 Alternative name(s): BRI1-kinase domain-interacting protein 105 Short name= BIP105
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.