



Recombinant Human Protein LDOC1L (LDOC1L)

Product Code	CSB-BP760786HU
Abbreviation	LDOC1L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ICC9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVQPQTSKAE SPALAASPNA QMDDVIDTLT SLRLTNSALR REASTLRAEK ANLTNMLESV MAELTLLRTR ARIPGALQIT PPISSITSNG TRPMTTPPTS LPEPFGSDPG RLAGFLMQMD RFMIFQASRF PGEAERVAFL VSRLTGAEAK WAIPHMQPDS PLRNYYQGFL AELRRTYKSP LRHARRAQIR KTSASNRAVR ERQMLCRQLA SAGTGPCPVH PASNGTSPAP ALPARARNL
Source	Baculovirus
Target Names	RTL6
Protein Names	Recommended name: Protein LDOC1L Alternative name(s): Leucine zipper protein down-regulated in cancer cells-like Mammalian retrotransposon-derived protein 6
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.