



Recombinant Human Protein LDOC1L (LDOC1L)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP760786HU |
| Abbreviation | LDOC1L |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6ICC9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MVQPQTSKAE SPALAASPNA QMDDVIDTLT SLRLTNSALR REASTLRAEK ANLTNMLESV MAELTLRTR ARIPGALQIT PPISSITSNG TRPMTTPPTS LPEPFSGDPG RLAGFLMQMD RFMIFQASRF PGEAERVAFL VSRLTGAEAK WAIPHMQPDS PLRNYYQGFL AELRRTYKSP LRHARRAQIR KTSASNRAVR ERQMLCRQLA SAGTGPCPVH PASNGTSPAP ALPARARNL |
| Source | Mammalian cell |
| Target Names | RTL6 |
| Protein Names | Recommended name: Protein LDOC1L Alternative name(s): Leucine zipper protein down-regulated in cancer cells-like Mammalian retrotransposon-derived protein 6 |
| Expression Region | 1-239 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |