



Recombinant *Oryza sativa* subsp. japonica Putative B3 domain-containing protein Os02g0455900 (Os02g0455900, LOC_Os02g25830)

Product Code	CSB-YP760842OFG
Abbreviation	Os02g0455900, LOC_Os02g25830
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6K3B1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAASPPLPTS IDGGQVLDDM EVVEMKYLFG KVLMPSDVSW DTEQLVIPDE HVGKLLDMVV MNRPEGGFFV VVVEDGEVTG KLWLFYWKR DDVHCLTKGW GCYAREKGLR AGDVSFFHS TACGRFFICC RCTCMSFLSL PTTSHRIHGS SVLPQPRAAQ EAHHPFSGHA TLCLGNKASD HSAPARHATA SLGCAAQPP QVPPTPRR RRRSMMVHPE PPEHTTDGMP VILESMALVS TPPVAKRVRVRL FGVYIDVPPL RPGGEATQDF NP
Source	Yeast
Target Names	Os02g0455900
Protein Names	Recommended name: Putative B3 domain-containing protein Os02g0455900
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.