



Recombinant Arabidopsis thaliana 14.7 kDa heat shock protein (HSP14.7)

Product Code	CSB-EP761203DOA
Abbreviation	HSP14.7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NLV0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSRNMEVNAG SSGEIPSPIR NRFQKSGSQA VYEVTTETKKS CVTRVDMPGC PESDLTYWVD ANNVHFFADE PAMPEYENAG RKYGGSMIFN PEAYDVKKTK VKLINGVLWI TVPKIPGKNA SINVKERILH Y
Source	E.coli
Target Names	HSP14.7
Protein Names	Recommended name: 14.7 kDa heat shock protein Short name= AtHsp14.7
Expression Region	1-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.