



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 58 (At5g18550)

Product Code	CSB-BP761207DOA
Abbreviation	At5g18550
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NPN3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MERYGGAGED ESRSDPSHEW SAQGTETGIE ASMWRLGLRG GGGGGETFPE RPDEPDCIYY LRTGVCGYGS RCRFNHPRNR APVLGGLRTE AGEFPERMGQ PVCQHFMRTG TCKFGASCKY HHPRQGGGGD SVTPVSLNYM GFPLRPGEKE CSYFMRTGQC KFGSTCRYHH PVPPGVQAPS QQQQQQLSAG PTMYPQLSQ TVPSSQQYGV VLARPQLLP SYVQSPYGYG QMVLPPGMVP YSGWNPYQAS VSAMPSPGTQ PSMGTSVYV ITPLSPSAPA YQSGPSSTGV SNKEQTFPQR PEQPECQYFM RTGDCKFGTS CRFHHPMEAA SPEASTLSHI GLPLRPGAVP CTHFAQHGIC KFGPACKFDH SLGSSSLSYS PSPSSLTDMV VAPYPSSLGT LAPSSSSDQC TELISSSSIE PITTTTGGSE TVAAGVSSMT SDVSHPEPAE TNKGDSASNE AKTSS
Source	Baculovirus
Target Names	At5g18550
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 58 Short name= AtC3H58 Alternative name(s): Zinc finger CCCH domain-containing protein ZFN-like 3
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.