



# Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 58 (At5g18550)

<b>Product Code</b>	CSB-EP761207DOA
<b>Abbreviation</b>	At5g18550
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NPN3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MERYGGAGED ESRSDPSHEW SAQGTETGIE ASMWRLGLRG GGGGGETFPE RPDEPDCIYY LRTGVCGYGS RCRFNHPRNR APVLGGLRTE AGEFPERMGQ PVCQHFMRTG TCKFGASCKY HHPRQGGGGD SVTPVSLNYM GFPLRPGEKE CSYFMRTGQC KFGSTCRYHH PVPPGVQAPS QQQQQQLSAG PTMYPQLQSQ TVPSSQQYGV VLARPQLLP SYVQSPYGYG QMVLPPGMVP YSGWNPYQAS VSAMPSPGTQ PSMGTSVYVYG ITPLSPSAPA YQSGPSSTGV SNKEQTFPQR PEQPECQYFM RTGDCKFGTS CRFHHPMEAA SPEASTLSHI GLPLRPGAVP CTHFAQHGIC KFGPACKFDH SLGSSSLSYS PSPSSLTDMV VAPYPSSLGT LAPSSSSDQC TELISSSSIE PITTTTGGSE TVAAGVSSMT SDVSHPEPAE TNKGDSASNE AKTSS
<b>Source</b>	E.coli
<b>Target Names</b>	At5g18550
<b>Protein Names</b>	Recommended name: Zinc finger CCCH domain-containing protein 58 Short name= AtC3H58 Alternative name(s): Zinc finger CCCH domain-containing protein ZFN-like 3
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.