



Recombinant Danio rerio Arginine/serine-rich coiled-coil protein 2 (rsrc2)

Product Code	CSB-MP761246DIL
Abbreviation	rsrc2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NWI1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAANDAE LHD S FSGKS QPGE NRKASRSSKH CSRSRSRSTE RKRKSGDKRH KRSHRSKEA RKKDSEKALK CQSGSEERLE FSDKGRDRLS EDTEERHRRK DKKTSRPRSH SRSRSKEKRH HNRNWDKRRS RSRSRSRDKK RRARSRSNSR SKHRHRSRSR SKSREKKKRI EKSRRKSRSP SISPVTFRGR NTAMDAQEAL ARRLERAKKL QEQKEKDMLE KWQHQEKAAA STPQC DPAAA APSPALNVAA LLASGTQVTP QIAMAAQMAA LQAKTLAETG IAVPSYYNPS AVNPMKFAEQ EKRRKMLWQG KKEGDNKSQT AELWEKLNFG NKDQNVKFRK LMGIKGEEEA ASSA AVNEDG LKTLQQQEEM FRNL DVQYEM ARSQTH TQRG MGLGFSSSFS SRGMDAV
Source	Mammalian cell
Target Names	rsrc2
Protein Names	Recommended name: Arginine/serine-rich coiled-coil protein 2
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.