



# Recombinant Human Interferon-inducible GTPase 5 (IRGC)

<b>Product Code</b>	CSB-EP761252HU-B
<b>Abbreviation</b>	IRGC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NXR0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATSKLPVVP GEEENTILMA KERLEALRTA FESGDLPQAA SHLQELLAST ESIRLEVGVT GESGAGKSSL INALRGLEAE DPGAALTGVM ETTMQPSPYP HPQFPDVTW DLPGAGSPGC PADKYLKQVD FSRDYFFLLV SPRRCGAVET RLAAEILCQG KKFYFVRTKV DEDLAATRTQ RPSGFREAAV LQEIRDHCAE RLREAGVADP RIFLVSNLSP ARYDFPTLVS TWEHDLPSHR RHAGLLSLPD ISLEALQKKK AMLQEQLVKT ALVLGVIQAL PVPGLAAAYD DALLIHSLRG YHRSFGLDDD SLAKLAEQVG KQAGDLRSVI RSPLANEVSP ETVLRLYSQS SDGAMRVARA FERGIPVFGT LVAGGISFGA VYTMLQGCLN EMAEDAQRVR IKALEDDEPQ PEVSLEVASD NGVEKGGSGE GGEEEAPLST CRKLGLLLKY ILDSWKKHDS EEK
<b>Source</b>	E.coli
<b>Target Names</b>	IRGC
<b>Protein Names</b>	Recommended name: Interferon-inducible GTPase 5 EC= 3.6.5.- Alternative name(s): Immunity-related GTPase cinema 1
<b>Expression Region</b>	1-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.