



Recombinant Human Interferon-inducible GTPase 5 (IRGC)

Product Code	CSB-YP761252HU
Abbreviation	IRGC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NXR0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MATSKLPVVP GEEENTILMA KERLEALRTA FESGDLPQAA SHLQELLAST ESIRLEVGVT GESGAGKSSL INALRGLEAE DPGAALTGVM ETTMQPSPYP HPQFPDVTW DLPGAGSPGC PADKYLKQVD FSRDYFFLLV SPRRCGAVET RLAAEILCQG KKFYFVRTKV DEDLAATRTQ RPSGFREAAV LQEIRDHCAE RLREAGVADP RIFLVSNLSP ARYDFPTLVS TWEHDLPSHR RHAGLLSLPD ISLEALQKKK AMLQEQLKLT ALVLGVIQAL PVPGLAAAYD DALLIHSLRG YHRSFGLDDD SLAKLAEQVG KQAGDLRSVI RSPLANEVSP ETVLRLYSQS SDGAMRVARA FERGIPVFGT LVAGGISFGA VYTMLQGCLN EMAEDAQRVR IKALEDDEPQ PEVSLEVASD NGVEKGGSGE GGEEAPLST CRKLGLLLKY ILDSWKKHDS EEK
Source	Yeast
Target Names	IRGC
Protein Names	Recommended name: Interferon-inducible GTPase 5 EC= 3.6.5.- Alternative name(s): Immunity-related GTPase cinema 1
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.