



Recombinant Danio rerio Protein MAK16 homolog (mak16)

Product Code	CSB-YP761255DIL
Abbreviation	mak16
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NYD4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MQHDDVIWDL IGNKSFCSYK VKTKTQQFCR NEYNITGLCN RSSCPLANSQ YATIREEKGG CFLYMKVIER AAFPSRMWEK VKLDRNYAKA LEQIDENLIY WPRFIRHKCK QRLTKITQYL IRIKLTTLKR QRKLVPLSRK VERREKRREE KALIAAQLN AIEKELLDRL KQGTYGDIYN FPINAFDKAM EKQDEESESE EEEEEEDEEE SGTREFIAED EFEESDLSDF EELNNLKDSS DEEVNDEEES SEEEEESEEE EEVKS KAKSK GKAPLKGPSR KKRAYVEIEY EQETEPAQKS KAT
Source	Yeast
Target Names	mak16
Protein Names	Recommended name: Protein MAK16 homolog Alternative name(s): MAK16-like protein Protein RBM13
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.