



Recombinant Danio rerio Exosome complex component MTR3 (exosc6)

Product Code	CSB-MP761268DIL
Abbreviation	exosc6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P0I8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MPVDTKRIRG PEESQSPLLF LSPDKAPKAP SSRQGVRGNG DVRPVFARCG LVSQAKGSAY IEAGNTKIIC SVYGPKETER RDETDMKTGR LVCDFRLAPF SCVKRGAWIQ GSEERDLSAT LMESLRPGVC LHRYPRESQID VNVMVLENDG SVLAHAVTCA SMALADAGIE MYDIVLGCTL RQSGNACLVD PSYAEECGSW QEGYGDNQGC VTLALLPNLN QVSGLNADGE MREDTLTEAM RTCMDGCHKL YPVVQQALTR AVKKKAPPPE K
Source	Mammalian cell
Target Names	exosc6
Protein Names	Recommended name: Exosome complex component MTR3 Alternative name(s): Exosome component 6 mRNA transport regulator 3 homolog
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.