



Recombinant Danio rerio 60S ribosomal protein L8 (rpl8)

Product Code	CSB-MP761270DIL
Abbreviation	rpl8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P0V6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	GRVIRGQQRK GAGSVFKAHV KHRKGAALKR HIDFAERHGY IKGIVKDIH DPGRGAPLAK VMFRDPYRFK KRTELFIAAE GIHTGQFIYC GKKAQLNIGN VLPVGTMPPEG TIVCCLEEKP GDRGKLARAS GNYATVISHN PETKKSrvKL PSGSKKVISS ANRAVVGvVA GGGRIDKPIL KAGRAYHkyK AKRNCWPRVR GVAMNPVEHP FGGGNHQHIG KPSTIRRDAP AGRKVGLIAA RRTGRLRGTK TVQDKEN
Source	Mammalian cell
Target Names	rpl8
Protein Names	Recommended name: 60S ribosomal protein L8
Expression Region	2-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.