



Recombinant Human PiggyBac transposable element-derived protein 2 (PGBD2)

Product Code	CSB-EP761285HU
Abbreviation	PGBD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P3X8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MASTSRDVIA GRGIHSKVKS AKLLEVLNAM EEEESNNNRE EIFIAPPDNA AGEFTDEDSG DEDSQRGAHL PGSVLHASVL CEDSGTGEDN DDLELQPAKK RQKAVVKPQR IWTKRDIRPD FGSWTASDPH IEDLKSQELS PVGLFELFFD EGTINFIVNE TNRYAWQKNV NLSLTAQELK CVLGILILSG YISYPRRRMF WETSPDSHHH LVADAIRRDR FELIFSYLHF ADNELDASD RFAKVRPLII RMNCNFQKHA PLEEFYSFGE SMCEYFGHRG SKQLHRGKPV RLGKIWCMT TSRGYLVWFE PSQGLTFTKP DRSLDLGGSM VIKFVDALQE RGFLPYHIFF DKVFTSVKLM SILRKKGVKA TGTVREYRTE RCPLKDPKEL KKMKRGSFDY KVDESEIIV CRWHDSSVVN ICSNAVGIEP VRLTSRHSGA AKTRTQVHQP SLVKLYQEKV GGVGRMDQNI AKYKVKIRGM KWYSSFIGYV IDAALNNAWQ LHRICCCDAQ VDLLAFRRYI ACVYLESNAD TTSQGRRSRR LETESRFDMI GHWIIHQDKR TRCALCHSQT NTRCEKCQKG VHAKCFREYH IR</p>
Source	E.coli
Target Names	PGBD2
Protein Names	Recommended name: PiggyBac transposable element-derived protein 2
Expression Region	1-592
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The piggyBac family of proteins, found in diverse animals, are transposases related to the transposase of the canonical piggyBac transposon from the moth, <i>Trichoplusia ni</i> . This family also includes genes in several genomes, including human, that appear to have been derived from the piggyBac transposons. This gene belongs to the subfamily of piggyBac transposable element derived (PGBD) genes. The PGBD proteins appear to be novel, with no obvious



relationship to other transposases, or other known protein families. The exact function of this gene is not known. Two transcript variants encoding different isoforms have been found for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.