



# Recombinant Danio rerio UPF0533 protein C5orf44 homolog (zgc:73187)

<b>Product Code</b>	CSB-MP761324DIL
<b>Abbreviation</b>	zgc:73187
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6PBY7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDVNQAKQEH LLALKVMRLT KPTLFTNMPV TCEDRDLPGD LFLRLMKDDP STVKGAETLI LGEMLTLPQN FGNIFLGETF SSYISVHNDS SQVVKDILVK ADLQTSSQRL NLSASNSAVS ELKPECCIDD VIHHEVKEIG THILVCAVSY TTQTGEKLYF RKFFKFQVLK PLDVKTKFYN AETDEVFLEA QIQNITTSPM FMEKVSLEPS MMYNVTELENN VASGDESSES TFGKMSYLQP LDTRQYLYCL KPKPEFAEKA GVIKGVTVIG KLDIVWKTNL GERGRLQTSQ LQRMAPGYGD VRLSLEFIPD TVDLEEPFDI TCKITNCSE TMDLLEMCN TRSVHWCGVS GRQLGKLSPTS ASLSIPLKLL SSVQGLQSSIS GLRLTDTFLK RTYEYDDIAQ VCVVCPYTSP ES
<b>Source</b>	Mammalian cell
<b>Target Names</b>	trappc13
<b>Protein Names</b>	Recommended name: UPF0533 protein C5orf44 homolog
<b>Expression Region</b>	1-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.