



Recombinant Human Zinc finger CCHC domain-containing protein 12 (ZCCHC12)

Product Code	CSB-BP761340HU
Abbreviation	ZCCHC12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PEW1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASIIARVGN SRRLNAPLPP WAHSMLRSLG RSLGPIMASM ADRNMKLFSG RVVPAQGEET FENWLTQVNG VLPDWNMSEE EKLKRLMKTL RGPAREVMRV LQATNPNSLV ADFLRAMKLV FGESESSVTA HGKFFNTLQA QGEKASLYVI RLEVQLQNAI QAGIIAEKDA NRTRLQQLLL GGELSRDLRL RLKDFLRMYA NEQERLPNFL ELIRMVREEE DWDDAFIKRK RPKRSESMVE RAVSPVAFQG SPPIVIGSAD CNVIEIDDTL DDSDEDVILV ESQDPPLPSW GAPPLRDRAR PQDEVVIDS PHNSRAQFPS TSGGSGYKNN GPGEMRRARK RKHTIRCSYC GEEGHSKETC DNESDKAQVF ENLIITLQEL THTEMERSRV APGEYNDFSE PL
Source	Baculovirus
Target Names	ZCCHC12
Protein Names	Recommended name: Zinc finger CCHC domain-containing protein 12 Alternative name(s): Smad-interacting zinc finger protein 1
Expression Region	1-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.