



# Recombinant Human Zinc finger CCHC domain-containing protein 12 (ZCCHC12)

<b>Product Code</b>	CSB-MP761340HU
<b>Abbreviation</b>	ZCCHC12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6PEW1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASIIARVGN SRRLNAPLPP WAHSMLRSLG RSLGPIMASM ADRNMKLFSG RVVPAQGEET FENWLTQVNG VLPDWNMSEE EKLKRLMKTL RGPAREVMRV LQATNPNSLV ADFLRAMKLV FGESESSVTA HGKFFNTLQA QGEKASLYVI RLEVQLQNAI QAGIIAEKDA NRTRLQQLLL GGELSRDLRL RLKDFLRMYA NEQERLPNFL ELIRMVREEE DWDDAFIKRK RPKRSESMVE RAVSPVAFQG SPPIVIGSAD CNVIEIDDTL DDSDEDVILV ESQDPPLPSW GAPPLRDRAR PQDEVVIDS PHNSRAQFPS TSGGSGYKNN GPGEMRRARK RKHTIRCSYC GEEGHSKETC DNESDKAQVF ENLIITLQEL THTEMERSRV APGEYNDFSE PL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZCCHC12
<b>Protein Names</b>	Recommended name: Zinc finger CCHC domain-containing protein 12 Alternative name(s): Smad-interacting zinc finger protein 1
<b>Expression Region</b>	1-402
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.