



Recombinant Ostreid herpesvirus 1

Uncharacterized protein ORF110 (ORF110)

Product Code	CSB-EP761402OBAD-B
Abbreviation	ORF110
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6R7C0
Product Type	Recombinant Protein
Immunogen Species	Ostreid herpesvirus 1 (isolate France) (OsHV-1) (Pacific oyster herpesvirus)
Purity	>85% (SDS-PAGE)
Sequence	MASKVYIDAF EDEWVLESDI QSTTTVTSDD KDEEMYYIID RDTRKPVGFK IPKDDEKKLF LFDDKFRYIF GAASIDYAYY SPLCKYKVIT SRTSPYLSHS VQYSRKTGRP GGREVCSSKEE IERIDSLFTV IDVNAMESKL DKVKDIREST IAEEYQIVLG KLKKDYIDFI MTGSELENRF ETCVKVKFLL EKSQKISRHE PRSFYSMVDR IAPHTGLIMQ NCNEDFSFKK AVDIAFFVSK QGSVLLNSHT KIATALKMIT N
Source	E.coli
Target Names	ORF110
Protein Names	Recommended name: Uncharacterized protein ORF110
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.