



Recombinant Human Dehydrogenase/reductase SDR family member 13 (DHRS13)

Product Code	CSB-BP761493HU
Abbreviation	DHRS13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UX07
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	PPCGG MGNLRGRTAV VTGANSIGIGK MTALELARRG ARVVLACRSQ ERGEAAAFDL RQESGNNEVI FMALDLASLA SVRAFATAFL SSEPRLDILI HNAGISSCGR TREFNLLLR VNHIGPFLLT HLLLPCLKAC APSRVVVVAS AAHCRGRLDF KRLDRPVVGW RQELRAYADT KLANVLFARE LANQLEATGV TCYAAHPGPV NSELFLRHVP GWLRPLLRLPL AWLVLRAPRG GAQTPLYCAL QEGIEPLSGR YFANCHVEEV PPAARDRAA HRLWEASKRL AGLGPGEDAE PDEDPQSEDS EAPSSLSTPH PEEPTVSQPY PSPQSSPDLS KMTHRIQAKV EPEIQLS
Source	Baculovirus
Target Names	DHRS13
Protein Names	Recommended name: Dehydrogenase/reductase SDR family member 13 EC=1.1.-.-
Expression Region	26-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.