



Recombinant Human Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2 (BAIAP2L2)

Product Code	CSB-EP761505HU-B
Abbreviation	BAIAP2L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UXY1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MAPEMDQFYR STMAIYKSIM EQFNPALENL VYLGNNYLRA FHALSEAAEV YFSAIQKIGE RALQSPTSQI LGEILVQMSD TQRHLNSDLE VVVQTFHGGL LQHMEKNTKL DMQFIKDSRQ HYELEYRHRA ANLEKCMSEL WRMERKRDKN VREMKESVNR LHAQMQUAFVS ESQRAAELEE KRRYRFLAEK HLLLSNTFLQ FFGRARGMLQ NRVLLWKEQS EASRSPSRAH SPGLLGALG PPYPSGRLTP TCLDMPPRPL GEFSSPRSRH GSGSYGTEPD ARPASQLEPD RRSLPRTPSA SSLYSGSAQS SRSNSFGERP GGGGGARRVR ALVSHSEGAN HTLLRFSAGD VVEVLVPEAQ NGWLYGKLEG SSASGWFPEA YVKALEEGPV NPMTPTVTPMT SMTSMSPMTP MNPGNELPSR SYPLRGSHSL DLLLDRPGNS IAPSEYWDGQ SRSRTPSRVP SRAPSPAPP LPSSRRSSMG STAVATDVKK LMSSEQYPPQ ELFPRGTNPF ATVKLRPTIT NDRSAPLIR </p>
Source	E.coli
Target Names	BAIAP2L2
Protein Names	Recommended name: Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2 Short name= BAI1-associated protein 2-like protein 2 Alternative name(s): Planar intestinal- and kidney-specific BAR domain protein Short name
Expression Region	1-529
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.