



# Recombinant Human Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2 (BAIAP2L2)

<b>Product Code</b>	CSB-EP761505HU
<b>Abbreviation</b>	BAIAP2L2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6UXY1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAPEMDQFYR STMAIYKSIM EQFNPALENL VYLGNNYLRA FHALSEAAEV YFSAIQKIGE RALQSPTSQI LGEILVQMSD TQRHLNSDLE VVVQTFHGGL LQHMEKNTKL DMQFIKDSRQ HYELEYRHRA ANLEKCMSEL WRMERKRDKN VREMKESVNR LHAQMQUAFVS ESQRAAELEE KRRYRFLAEK HLLLSNTFLQ FFGRARGMLQ NRVLLWKEQS EASRSPSRAH SPGLLGPALG PPYPSGRLTP TCLDMPPRPL GEFSSPRSRH GSGSYGTEPD ARPASQLEPD RRSLPRTPSA SSLYSGSAQS SRSNSFGERP GGGGGARRVR ALVSHSEGAN HTLLRFSAGD VVEVLVPEAQ NGWLYGKLEG SSASGWFPEA YVKALEEGPV NPMTPTVTPMT SMTSMSPMTP MNPGNELPSR SYPLRGSHSL DDLLDRPGNS IAPSEYWDGQ SRSRTPSRVP SRAPSPAPP LPSSRRSSMG STAVATDVKK LMSSEQYPPQ ELFPRGTNPF ATVKLRPTIT NDRSAPLIR
<b>Source</b>	E.coli
<b>Target Names</b>	BAIAP2L2
<b>Protein Names</b>	Recommended name: Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2 Short name= BAI1-associated protein 2-like protein 2 Alternative name(s): Planar intestinal- and kidney-specific BAR domain protein Short name
<b>Expression Region</b>	1-529
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.