



Recombinant Human Transcription factor AP-2-epsilon (TFAP2E)

Product Code	CSB-EP761516HU-B
Abbreviation	TFAP2E
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6VUC0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLVHTYSAME RPDGLGAAAG GARLSSLPQA AYGPAPPLCH TPAATAAAEF QPPYFPPYP QPPLPYGQAP DAAAAPHLA GDPYGGGLAPL AQPQPQA AAPRAAARAH EEPGLLAPP ARALGLDPRR DYATAVPRLL HGLADGAHGL ADAPLGLPGL AAAPGLEDLQ AMDEPGMSLL DQSVIKKVPI PSKASSLSAL SLAKDSLGG ITNPGEVFCV VPGRLSLLSS TSKYKVTVE VQRRLSPEEC LNASLLGGVL RRAKSKNGGR CLRERLEKIG LNLPAARRKA ANVTLLTSLV EGEAVHLARD FGYVCETFP AKAAAEYLCR QHADPGELHS RKSMLLAQK ICKEFADLMA QDRSPLGNSR PALILEPGVQ SCLTHFSLIT HGFGGPAICA ALTAFQNYLL ESLKGLDKMF LSSVGSNGHGE TKASEKDAKH RK
Source	E.coli
Target Names	TFAP2E
Protein Names	Recommended name: Transcription factor AP-2-epsilon Short name= AP2-epsilon Alternative name(s): Activating enhancer-binding protein 2-epsilon
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.