



# Recombinant Human Chordin-like protein 2 (CHRDL2)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP761533HU   |
| <b>Abbreviation</b>      | CHRDL2   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q6WN34   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | RARPD MFCLFHGKRY SPGESWHPYL EPQGLMYCLR CTCSEGAHVS<br>CYRLHCPPVH CPQPVTPEQQ CCPKCVEPHT PSGLRAPPKS<br>CQHNGTMYQH GEIFSAHELF PSRLPNQCVL CSCTEGQIYC GLTTCPEPGC<br>PAPLPLPDSC CQACKDEASE QSDEEDSVQS LHGVRHPQDP<br>CSSDAGRKRK PGTPAPTGLS APLSFIPRHF RPKGAGSTTV KIVLKEKHKK<br>ACVHGGKTYS HGEVWHPAFR AFGPLPCILC TCEDGRQDCQ RVTCPTEYPC<br>RHPEKVAGKC CKICPEDKAD PGHSEISSTR CPKAPGRVLV HTSVSPSPDN<br>LRRFALEHEA SDLVEIYLWK LVKGFHFLTQ IKKVRKQDFQ KEAQHFRLLA<br>GPHEGHWNVF LAQTLELKVT ASPDKVTKT |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | CHRDL2   |
| <b>Protein Names</b>     | Recommended name: Chordin-like protein 2 Alternative name(s): Breast tumor novel factor 1 Short name= BNF-1 Chordin-related protein 2  |
| <b>Expression Region</b> | 26-429   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.