



Recombinant Human Hydroxysteroid dehydrogenase-like protein 2 (HSDL2)

Product Code	CSB-EP761576HU
Abbreviation	HSDL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YN16
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLPNTGRLAG CTVFITGASR GIGKAIALKA AKDGANIVIA AKTAQPHPKL LGTIYTAAEE IEAVGGKALP CIVDVRDEQQ ISAAVEKAIK KFGGIDILVN NASAI SLTNT LDTPTKRLDL MMNVNTRGTY LASKACIPYL KSKVAHILN ISPPLNLPV WFKQH CAYTI AKYGMSMYVL GMAEEFKGEI AVNALWPKTA IHTAAMDMLG GPGIESQCRK VDIADAAYS IFQKPKSFTG NFVIDENILK EEGIENFDVY AIKPGHPLQP DFFLDEYPEA VSKKVESTGA VPEFKEEKLQ LQPKPRGAV EETFRIVKDS LSDDVVKATQ AIYLFELSGE DGGTWFLDLK SKGGNVGYGE PSDQADVMS MTTDDFVKMF SGKLPKPTMAF MSGKLIKIKGN MALAIKLEKL MNQMNARL
Source	E.coli
Target Names	HSDL2
Protein Names	Recommended name: Hydroxysteroid dehydrogenase-like protein 2 EC= 1.-.-.-
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.