



Recombinant *Oryza sativa* subsp. japonica Putative B3 domain-containing protein Os08g0333500 (Os08g0333500, LOC_Os08g24470)

Product Code	CSB-EP761593OFG-B
Abbreviation	Os08g0333500, LOC_Os08g24470
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YTQ4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTVELEKIAG SFFISKGWKT FVHRTGLLSG QYIRFQVLTP SKINVLLFDK KKDSKLPMP SSKKQIKTAP KRSTGITIND MPTSKHASML ISHTSNKETS SDSRTESMTD IPSSSDNSGE TTRSFDDLFC CARNTAVTPD IKNYISIIGQ FLQRSSKFYI VTMNNTFMKQ DRLCKLSLLT NMLMCDVLTD TCRFTLPSFL MAFFQNTVQT QTLTTCADSP YPYEHLRKTE PACREIDEVT LTMCAHSPYP YEYLRKTELA CREIDEVTTG ASLSTNILPP YTSSPLEAIS GGEDIWWLAA SRMGWRGDVA AVQADSV DGN NSPFLCARSP TIGNGGADQP YLMVGKPHGW PCWGWDGAEM WLWCQSKNSG HCHIGWPVGV LVKGGRIYPG GGLPSKELAG MIGPYASGDW FELKRELPK A
Source	E.coli
Target Names	Os08g0333500
Protein Names	Recommended name: Putative B3 domain-containing protein Os08g0333500
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.