



# Recombinant *Oryza sativa* subsp. japonica Putative B3 domain-containing protein Os08g0333500 (Os08g0333500, LOC\_Os08g24470)

<b>Product Code</b>	CSB-EP761593OFG
<b>Abbreviation</b>	Os08g0333500, LOC_Os08g24470
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6YTQ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTVELEKIAG SFFISKGWKT FVHRTGLLSG QYIRFQVLTP SKINVLLFDK KKDSKLPMP SSKKQIKTAP KRSTGITIND MPTSKHASML ISHTSNKETS SDSRTESMTD IPSSSDNSGE TTRSFDDLFC CARNTAVTPD IKNYISIIGQ FLQRSSKFYI VTMNNTFMKQ DRLCKLSLLT NMLMCDVLTD TCRFTLPSFL MAFFQNTVQT QTLTTCADSP YPYEHLRKTE PACREIDEVT LTMCAHSPYP YEYLRKTELA CREIDEVTTG ASLSTNILPP YTSSPLEAIS GGEDIWWLAA SRMGWRGDVA AVQADSV DGN NSPFLCARSP TIGNGGADQP YLMVGKPHGW PCWGWDGAEM WLWCQSKNSG HCHIGWPVGV LVKGGRIYPG GGLPSKELAG MIGPYASGDW FELKRELPK A
<b>Source</b>	E.coli
<b>Target Names</b>	Os08g0333500
<b>Protein Names</b>	Recommended name: Putative B3 domain-containing protein Os08g0333500
<b>Expression Region</b>	1-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.