



Recombinant *Oryza sativa* subsp. *japonica* COBRA-like protein 6 (BC1L7)

Product Code	CSB-MP761611OFG
Abbreviation	BC1L7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z4G7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	YDPLDPNGN ITIKWDVMSW TPDGYVAMVT INNYQTYRQI MAPGWTVGWT WARQEVISM VGAQATDQGD CSRFKANLPH CCRRTPAVVD LLPGVPYNQQ IANCCRGGVL PAYGQAPSAA AAFQVSVGQ AGTTNRTVRL PRNFTLLGPG PGYTCGRARV VPSTVFLTAD RRRKTQALMT WNVCTYSQH LASKYPSCCV SFSSFYNDTI VPCAKCACGC DAHKPCVRSE RDGKRLAVTG KKHDANANAH GRGNGVAAAA MAAPLLQCTT HMCPVRVHWH VKLNYREYWR AKITIVNFNY RMNYTGWTLV AQHPNLDNIT EVFSFDYKPV VSYGSINDTA MFYGLKYFND QLMEAGPHGN VQSEVLMRKD ARTFTFRQGW AFPRKVYFNG DECQMPPPS YPYLPN
Source	Mammalian cell
Target Names	BC1L7
Protein Names	Recommended name: COBRA-like protein 6 Alternative name(s): Protein BRITTLE CULM1-like 7
Expression Region	22-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.