



Recombinant *Oryza sativa* subsp. japonica Probable 6-phosphogluconolactonase 3, chloroplastic (Os08g0547100, LOC_Os08g43370)

Product Code	CSB-MP761615OFG
Abbreviation	Os08g0547100, LOC_Os08g43370
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z9C3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	AAARKKLLI FDAEEYLAES LAKYTARLSG EAVAERGAFT VALSGGSLIK ALRKLTESPY LEAVEWSKWH VFVVDERVVP KDHADSNYKL AMDGLLSKVP IPASQIYAIN DTLSAEGAAD EYETCLKQLV NDGVVAISEV TGFPKLDLML LGMGPDGHVA SLFPGHPVVN ENLKWVSYIK DSPKPPPERI TFTFPLVNSS AHIALVVTGA GKAGAVHKAF SDKQSSSDLL PVEMVSQQEG VLTWFTDKPA VSMLSSI
Source	Mammalian cell
Target Names	Os08g0547100
Protein Names	Recommended name: Probable 6-phosphogluconolactonase 3, chloroplastic Short name= 6PGL 3 EC= 3.1.1.31
Expression Region	72-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.