



Recombinant Human ADAMTS-like protein 5 (ADAMTSL5)

Product Code	CSB-BP761624HU
Abbreviation	ADAMTSL5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZMM2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QGPGEWTP WVSWTRCSSS CGRGVSVRSR RCLRLPGEEP CWGDSHEYRL CQLPDCPPGA VPFRLDQCAL YNGRPVLGTQ KTYQWVPFHG APNQCDLNCL AEGHAFYHSF GRVLDGTACS PGAQGVCVAG RCLSAGCDGL LGSGALEDRC GRCGGANDSC LFVQRVFRDA GAFAGYWNVT LIPEGARHIR VEHRSRNHLA LMGGDGRYVL NGHWVVSPPG TYEAAGTHVV YTRDTGPQET LQAAGPTSHD LLLQVLLQEP NPGIEFEFWL PRERYSPFQA RVQALGWPLR QPQPRGVEPQ PPAAPAVTPA QTPTLAPDPC PPCPDTRGRA HRLHYCGSD FVFQARVLGH HHQAQETRYE VRIQLVYKNR SPLRAREYVW APGHCPCPML APHRDYLMAY QRLVSPDGTQ DQLLLPHAGY ARPWSPAEDS RIRLTARRCP G
Source	Baculovirus
Target Names	ADAMTSL5
Protein Names	Recommended name: ADAMTS-like protein 5 Short name= ADAMTSL-5 Alternative name(s): Thrombospondin type-1 domain-containing protein 6
Expression Region	43-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.