



Recombinant Human Homeobox protein DBX2 (DBX2)

Product Code	CSB-BP761629HU
Abbreviation	DBX2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZNG2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLPSAVAAHA GAYWDVVASS ALLNLPAAPG FGNLGKSFLI ENLLRVGGAP TPRLQPPAPH DPATALATAG AQLRPLPASP VPLKLCPAAE QVSPAGAPYG TRWAFQVLSP SADSARLPGR APGDRDCTFQ PSAPAPSKPF LLSTPPFYSA CCGGSCRRA SSTAFFREES MLPLLTQDSN SKARRGILRR AVFSEDQRKA LEKMFQKQKY ISKTDRKKLA INLGLKESQV KIWFQNRMMK WRNSKEKEVL SNRCIQEVGL QEDPLRSAL GFPSPCPSIW DVPQQHSSPR WRENSPEPSE RLIQESSGAP PPEANSLQGA LYLCSEEEAG SKGVLTGAV
Source	Baculovirus
Target Names	DBX2
Protein Names	Recommended name: Homeobox protein DBX2 Alternative name(s): Developing brain homeobox protein 2
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.