



Recombinant Chicken Nuclear transcription factor Y subunit beta (NFYB)

Product Code	CSB-MP015774CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P25207
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MDGDSSTTDA SQLGIAGDYI GGSHYVIQPH DDTEDSMNDH EDTNGSKESF REQDIYLPYA NVARIMKNAI PQTGKIAKDA KECVQECVSE FISFITSEAS ERCHQEKRKT INGEDILFAM STLGFDYSYVE PLKLYLQKFR EAMKGEKGIG GTVTTGDGLS EELTEEAFTN QLPAGLITTD GQQQNVMVYT TSYQQISGVQ QIQFS
Source	Mammalian cell
Target Names	NFYB
Protein Names	Recommended name: Nuclear transcription factor Y subunit beta Alternative name(s): CAAT box DNA-binding protein subunit B Nuclear transcription factor Y subunit B Short name= NF-YB
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is one subunit of a trimeric complex, forming a highly conserved transcription factor that binds with high specificity to CCAAT motifs in the promoter regions in a variety of genes. This gene product, subunit B, forms a tight dimer with the C subunit, a prerequisite for subunit A association. The resulting trimer binds to DNA with high specificity and affinity. Subunits B and C each contain a histone-like motif. Observation of the histone nature of these subunits is supported by two types of evidence; protein sequence alignments and experiments with mutants.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.