



Recombinant *Ashbya gossypii* Serine/threonine-protein kinase BUD32 (BUD32)

Product Code	CSB-EP762168DOT-B
Abbreviation	BUD32
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q74Z75
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSAEVISRVS QIYGEGITLT PISQGAEAVI FTTTTHPYLP STEGSSKYVV KYRPPKTYRH PSIDVALTKQ RTLGEARLLG KLQIQGLRV PKLLACDVYN GSIWMEFLGE DLPESRGFSN LKNFLWMCAS EPYSPVEAT VRDVGCCQIGL LHWNDYCHGD LTTSNIVLVR SGAEWAPHLI DFGLGSNSNL VEDKGVLDLYV LERALASTHS AFAEKYFEWL MEGFTSVYQS HGEKGEQKLG ELLRRFEEVR LRGRKRSMLG
Source	<i>E.coli</i>
Target Names	BUD32
Protein Names	Recommended name: Serine/threonine-protein kinase BUD32 EC= 2.7.11.1
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.